



Product Information



Sterile hypodermic needle for single use

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect¹.

| Description | Raw materials |
|-------------|------------------------------|
| 1 Hub | Polypropylene + master batch |
| 2 Cannula | Stainless steel |
| Adhesive | Epoxy glue |
| Lubricant | Silicone oil |

Does not contain components made of PVC or natural rubber latex.



Sterilization

With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶.



Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the packaging
- Do not store at extreme temperature and humidity



Precautions

- Check integrity of the package before use
- Do not reuse
- Do not use if unit package is damaged
- Dispose of safely after single use to avoid risk of infection

¹ Data on file

Your needles of choice for a wide range of applications

Product range

| Product Code | Colour | Cannula ø | Cannula Length | Bevel |
|--------------|--------|-------------|----------------|-------------|
| AN*1838R1 | Pink | 18G 1.2 mm | 1 ½ inch 38 mm | Regular 11° |
| AN*1838S1 | Pink | 18G 1.2 mm | 1 ½ inch 38 mm | Short 17° |
| AN*1925R1 | Cream | 19G 1.1 mm | 1 inch 25 mm | Regular 11° |
| AN*1938R1 | Cream | 19G 1.1 mm | 1 ½ inch 38 mm | Regular 11° |
| AN*1938S1 | Cream | 19G 1.1 mm | 1 ½ inch 38 mm | Short 17° |
| AN*1950R1 | Cream | 19G 1.1 mm | 2 inch 50 mm | Regular 11° |
| AN*2025R1 | Yellow | 20G 0.9 mm | 1 inch 25 mm | Regular 11° |
| AN*2038R1 | Yellow | 20G 0.9 mm | 1 ½ inch 38 mm | Regular 11° |
| AN*2116R1 | Green | 21G 0.8 mm | ⅝ inch 16 mm | Regular 11° |
| AN*2125R1 | Green | 21G 0.8 mm | 1 inch 25 mm | Regular 11° |
| AN*2138R1 | Green | 21G 0.8 mm | 1 ½ inch 38 mm | Regular 11° |
| AN*2150R1 | Green | 21G 0.8 mm | 2 inch 50 mm | Regular 11° |
| AN*2232R1 | Black | 22G 0.7 mm | 1 ¼ inch 32 mm | Regular 11° |
| AN*2238R1 | Black | 22G 0.7 mm | 1 ½ inch 38 mm | Regular 11° |
| AN*2316R1 | Blue | 23G 0.6 mm | ⅝ inch 16 mm | Regular 11° |
| AN*2325R1 | Blue | 23G 0.6 mm | 1 inch 25 mm | Regular 11° |
| AN*2332R1 | Blue | 23G 0.6 mm | 1 ¼ inch 32 mm | Regular 11° |
| AN*2425R1 | Purple | 24G 0.55 mm | 1 inch 25 mm | Regular 11° |
| AN*2516R1 | Orange | 25G 0.5 mm | ⅝ inch 16 mm | Regular 11° |
| AN*2525R1 | Orange | 25G 0.5 mm | 1 inch 25 mm | Regular 11° |
| AN*2613R1 | Brown | 26G 0.45 mm | ½ inch 13 mm | Regular 11° |
| AN*2623R1 | Brown | 26G 0.45 mm | ⅞ inch 23 mm | Regular 11° |
| AN*2716R1 | Grey | 27G 0.4 mm | ⅝ inch 16 mm | Regular 11° |
| AN*2719R1 | Grey | 27G 0.4 mm | ¾ inch 19 mm | Regular 11° |
| AN*3013R1 | Yellow | 30G 0.3 mm | ½ inch 13 mm | Regular 11° |

Packaging

Single unit Blister pack (film/paper)

Unit box 100 pieces

Dimensions 105 x 84 x 86 mm

For item numbers AN*1950R1 and

AN*2150R1 dimensions 105 x 84 x 110 mm

Shipping carton 10.000 pieces

Dimensions 450 x 440 x 450 mm

For item numbers AN*1950R1 and

AN*2150R1 5.000 pieces,

dimensions 545 x 240 x 450 mm

Regulatory Status

Class IIa, sterile product following

MDD 93/42/EEC

Notified body TÜV Product Service

GmbH, TÜV SUD Group 0123

Manufacturer Zhejiang Kindly Medical

Devices Co. Ltd. (China)

EC representative Shanghai International

Corp. GmbH (Europe) - Hamburg, Germany

Distributor TERUMO EUROPE N.V.

Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a wide range from 18G up to 30G

The outside diameter of the needle can be recognized by the ISO colour-coded hub.



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Your Terumo partner